


**ET 6006 MINI - MATIC - 3 THREE MODULES + SYRINGE
ON MINI CART**

ET 6005 MINI - MATIC - 2 Two modules + Syringe
ET 6006 MINI - MATIC - 3 Three modules + Syringe


REF. NO.	DESCRIPTION	PART NO.	REF. NO.	DESCRIPTION
1	Air indicator light	UN 0300	37	Air return from toggle - master
2	Toggle valve , gray , master on - off	DC 7016	38	Air to regulator
3	Needle valve , water adjustment	UN 0200	39	Air to foot control
4	Toggle valve , gray , air on - off for water relay	DC 7016	40	Air to connector assembly
5	Needle valve , air flow adjustment	UN 0200	41	Air to automatic handpiece holder
6	Toggle valve , gray , air on - off spray	DC 7016		
7	Rectangular pressure gauge , 60 PSI	DC 7270	43	Air to syringe
8	Automatic handpiece holder No 1,gray	DC 4561		
9	Automatic handpiece holder No 2,gray	DC 4561	45	Air to control block from foot control
10	Automatic handpiece holder No 3,gray	DC 4561	46	Air to water relay command
11	Control block No 1	UN 0611	47	Air spray
12	Control block No 2	UN 0621	48	Air spray to block from chip blower
13	Control block No 3	UN 0631	49	Air to rectangular pressure gauge
14	Connector assembly	UN 0180	50	Water to water filter
15	Water relay	UN 0740	51	Water water relay
16	Needle valve , water adjustment for syringe	UN 0200	52	Water to block from water relay
			53	Water to syringe
			54	Air to toggle - master
19	Syringe	DC	55	Water adjustment
20	Foot control with chip blower	FC 2000	56	Drive pressure adjustment
			57	Air spray adjustment
			58	Water to turbine
			59	Air command to pick - up
24	Check valve 1/8 " O.D. tubing	UN 0940	60	Drive to turbine
25	End cap No 1	UN 0481	61	Air spray to turbine
26	Stud - bolt M4 , SST	UN 0501		
27	Pic - up	UN 0441		
28	Water filter	JB 0100		
29	Air filter miniature for regulator	WI F002		
30	Air gauge	WI G227		
31	Floor air valve	JB 0300		
32	Floor water valve	JB 0200		
33	Manual ball valve SA 1/2 "	SA 8000		
34	Air regulator 0 - 125 PSI	WI R001		
35	Adapter m/w , PUR	COA 12		
36	Adapter - borden	COA 36		

**ET 6006 MINI MATIC - 3, THREE MODULES + SYRINGE
WITH INTEGRATED CLEAN WATER SYSTEM
on MINI CART SCHEMA**



**ET 6006 MINI MATIC - 3 , THREE MODULES + SYRINGE
WITH INTEGRATED CLEAN WATER SYSTEM
on MINI CART SCHEMA
(Continued)**

REF. NO.	DESCRIPTION	PART NO.	REF. NO.	DESCRIPTION
1	Air indicator light	UN 0300		
2	Toggle valve , gray , master on - off	DC 7016	37	Air return from toggle - master
3	Needle valve , water adjustment	UN 0200	39	Air to foot control
4	Toggle valve , gray , air on - off for water relay	DC 7016	40	Air to connector assembly
5	Needle valve , air flow adjustment	UN 0200	41	Air to automatic handpiece holder
		DC 7016	43	Air to syringe
7	Rectangular pressure gauge , 60 PSI	DC 7270	45	Air to control block from foot control
8	Automatic handpiece holder No 1,gray	DC 5960	46	Air to water relay command
9	Automatic handpiece holder No 2,gray	DC 5960	47	Air spray
10	Automatic handpiece holder No 3,gray	DC 5960	48	Air spray to block from chip blower
11	Control block No 1	UN 0611	49	Air to rectangular pressure gauge
12	Control block No 2	UN 0621	50	Water to connector assembly
13	Control block No 3	UN 0631	51	Water to water relay
14	Connector assembly	UN 0180	52	Water to block from water relay
15	Water relay	UN 0740	53	Water to syringe
16	Needle valve , water adjustment for syringe	UN 0200	54	Air to toggle - master
17	Tee 1/8 " O.D. tubing	UN 0920	55	Water adjustment
19	Syringe	DC	56	Drive pressure adjustment
20	Foot control with chip blower	FC 2000	57	Air spray adjustment
24	Check valve 1/8 " O.D. tubing	UN 0940	58	Water to turbine
25	End cap No 1	UN 0481	59	Air command to pick - up
26	Stud - bolt M4 , SST	UN 0501	60	Drive to turbine
27	Pic - up	UN 0441	61	Air spray to turbine
29	Air filter miniature for regulator	WI F002		
30	Air gauge	WI G227		
31	Clean water bottle 1.5 liter	UNC 030		
32	Non - relieving regulator for air	DC 7200		
34	Air regulator 0 - 125 PSI	WI R001		
35	Adapter m/w , PUR	COA 12		
36	Adapter - borden	COA 36		



MINI - MATIC - 2 ET 6002 Two automatic modules
MINI - MATIC - 3 ET 6003 Three automatic modules
 + Syringe + Chip blower

**MINI - MATIC - 2 ET 6002 Two automatic modules
 MINI - MATIC - 3 ET 6003 Three automatic modules
 + Syringe + Chip blower**

REF. NO.	DESCRIPTION	PART NO.	REF. NO.	DESCRIPTION
1	Air indicator light	UN 0300	37	Air return from toggle - master
2	Toggle valve , gray , master on - off	DC 7016	38	Air to regulator
3	Needle valve , water adjustment	UN 0200	39	Air to foot control
4	Toggle valve , gray , air on - off for water relay	DC 7016	40	Air to connector assembly
5	Needle valve , air flow adjustment	UN 0200	41	Air to automatic handpiece holder
6	Toggle valve , gray , air on - off spray	DC 7016		
7	Rectangular pressure gauge , 60 PSI	DC 7270	43	Air to syringe
8	Automatic handpiece holder No 1,gray	DC 4561		
9	Automatic handpiece holder No 2,gray	DC 4561	45	Air to control block from foot control
10	Automatic handpiece holder No 3,gray	DC 4561	46	Air to water relay command
11	Control block No 1	UN 0611	47	Air spray
12	Control block No 2	UN 0621	48	Air spray to block from chip blower
13	Control block No 3	UN 0631	49	Air to rectangular pressure gauge
14	Connector assembly	UN 0180	50	Water to water filter
15	Water relay	UN 0740	51	Water water relay
16	Needle valve , water adjustment for syringe	UN 0200	52	Water to block from water relay
			53	Water to syringe
			54	Air to toggle - master
19	Syringe	DC	55	Water adjustment
20	Foot control with chip blower	FC 2000	56	Drive pressure adjustment
			57	Air spray adjustment
			58	Water to turbine
			59	Air command to pick - up
24	Check valve 1/8 " O.D. tubing	UN 0940	60	Drive to turbine
25	End cap No 1	UN 0481	61	Air spray to turbine
26	Stud - bolt M4 , SST	UN 0501		
27	Pic - up	UN 0441		
28	Water filter	JB 0100		
29	Air filter miniature for regulator	WI F002		
30	Air gauge	WI G227		
31	Floor air valve	JB 0300		
32	Floor water valve	JB 0200		
33	Manual ball valve SA 1/2 "	SA 8000		
34	Air regulator 0 - 125 PSI	WI R001		
35	Adapter m/w , PUR	COA 12		
36	Adapter - borden	COA 36		