


# ET 041S UNO Matic - BOX AUTOMATIC


## 1 HP. WITH WATER, 1 HP. DRY



REF. NO.	DESCRIPTION	PART NO.
1	CONTROL BLOCK	UN 0402
2	NEEDLE VALVE, WATER	UN 0200
3	TOGGLE VALVE, BLACK, AIR ON-OFF FOR WATER RELAY	DC 7015
4	NEEDLE VALVE, AIR	UN 0200
5	TEE 1/4 x 1/4 x 1/4	UN 0910
6	TEE 1/4 x 1/4 x 1/8	UN 0930
7	FOOT CONTROL	FC 1000
8	AIR MICROSWITCH N.O.	UN 0422
9	WATER RELAY	UN 0740
10	ADAPTER M/W, PUR - I	COA 12
11	ADAPTER M/W, PUR - II	COA 12
12	STRINGE	DCI
13	AIR REGULATOR 0 - 125 PSI	WI R001
14	AIR FILTER FOR REGULATOR	WI F002
15	AIR GAUGE	WI G227

1/4" WATER	
1/8" WATER	
1/4" AIR	
1/8" AIR	

# ET 042S UNO Matic - BOX AUTOMATIC 2 HP. WITH WATER




REF. NO.	DESCRIPTION	PART NO.
1	CONTROL VALVE, WATER	UN 0492
2	NEEDLE VALVE, WATER I	UN 0200
3	NEEDLE VALVE, WATER II	UN 0200
4	TOGGLE VALVE, BLACK, AIR ON-OFF FOR WATER RELAY I	DC 7015
5	TOGGLE VALVE, BLACK, AIR ON-OFF FOR WATER RELAY II	DC 7015
6	NEEDLE VALVE, AIR I	UN 0200
7	NEEDLE VALVE, AIR II	UN 0200
8	AIR MICROSWITCH N.O.	UN 0422
9	WATER RELAY	UN 0749
10	ADAPTER M/W, PUR - I	COA 12
11	ADAPTER M/W, PUR - H	COA 12
12	SYRINGE	DCI
13	FOOT CONTROL	FC 1006
14	TEE 1/4 x 1/4 x 1/4	UN 0910
15	TEE 1/8 x 1/4 x 1/8	UN 0930
16	TEE 1/8 x 1/8 x 1/8	UN 0926
17	AIR REGULATOR 0 - 125 PSI	WI R001
18	AIR FILTER FOR REGULATOR	WI F002
19	AIR GAUGE	WI G227

27


27

# ET 4400 UNO Matic-BOX AUTOMATIC 1 HP WITH WATER, 1 HP DRY WITH INTEGRATED CLEAN WATER SYSTEM



REF NO.	DESCRIPTION	PART NO.
1	CONTROL BLOCK	UN 0402
2	NEEDLE VALVE, WATER	UN 0200
3	TOGGLE VALVE, BLACK, AIR ON-OFF FOR WATER RELAY	DC 7015
4	NEEDLE VALVE, AIR	UN 0200
5	TEE 1/4 x 1/4 x 1/4	UN 0910
6	TEE 1/4 x 1/4 x 1/8	UN 0930
7	FOOT CONTROL	FC 1000
8	AIR MICROSWITCH N.O.	UN 0422
9	WATER RELAY	UN 0740
10	ADAPTER M/W, PUR - I	COA 12
11	ADAPTER M/W, PUR - II	COA 12
12	SYRINGE	
13	AIR REGULATOR 0-125 PSI	WT R001
14	AIR FILTER	WT F002
15	AIR GAUGE	WT G227
16	TOGGLE VALVE, AIR ON-OFF FOR NON-RELIEVING REGULATOR	DC 7015
17	NON-RELIEVING REGULATOR FOR AIR	DC 7200
18	CLEAN WATER BOTTLE 1.5 LITER	

# ET 4401 UNO MATIC - BOX AUTOMATIC 2 HP WITH WATER WITH INTEGRATED CLEAN WATER SYSTEM



R/F NO.	DESCRIPTION	PART NO.
1	CONTROL BLOCK	UN 0482
2	NEEDLE VALVE, WATER I	UN 0200
3	NEEDLE VALVE, WATER II	UN 0200
4	TOGGLE VALVE, BLACK, AIR ON-OFF FOR WATER RELAY I	DC 7015
5	TOGGLE VALVE, BLACK, AIR ON-OFF FOR WATER RELAY II	DC 7015
6	NEEDLE VALVE, AIR I	UN 0200
7	NEEDLE VALVE, AIR II	UN 0200
8	AIR MICROSWITCH N.O.	UN 0422
9	WATER RELAY	UN 0740
10	ADAPTER M/W, PUR - I	COA 12
11	ADAPTER M/W, PUR - II	COA 12
12	STRINCE	
13	FOOT CONTROL	PC 1000
14	TEE 1/4 X 1/4 X 1/4	UN 0910
15	TEE 1/4 X 1/4 X 1/8	UN 0930
16	TEE 1/8 X 1/8 X 1/8	UN 0920
17	AIR REGULATOR 0-125 PSI	WT R091
18	AIR FILTER FOR REGULATOR	WT F082
19	AIR GAUGE	WT G227
20	TOGGLE VALVE, AIR ON-OFF FOR NON-RELIEVING REGULATOR	DC 7015

21	NON-RELIEVING REGULATOR FOR AIR	DC 7200
22	CLEAN WATER BOTTLE 1.5 LITER	

1/8" WATER	
1/4" AIR	
1/8" AIR	