

# ETI

## DENTAL INDUSTRIES LTD.


New Tel. For Factory and show Room 972-3-5799760, 972-3-5799759

16, HAKISHON Bney-Brak 51203 ISRAEL

Mail Address P.O.BOX: 1400 Bney-Brak 51114 Israel.

FAX. No. 972-3-5799758

**ADAPTOR FOR CHAIR TECNODENT ECO**  
**WITH POST FOR UNIT, LIGHT, CUSPIDOR**



# ETI

## DENTAL INDUSTRIES LTD.

New Tel. For Factory and show Room 972-3-5799760, 972-3-5799759


16, HAKISHON Bney-Brak 51203 ISRAEL

Mall Address P.O.BOX: 1400 Bney-Brak 51114 Israel.

FAX. No. 972-3-5799758

### HEIGHT OF POST FOR UNIT & CUSPIDOR ON DIFFERENT CHAIR

REF. NO.	CHAIR	POST for UNIT (cm)
1	TECNODENT ECO	42
2	TECNODENT LINDA	37
3	JUMBO SDS	42
4	PELTON & CRANE	37



REF. NO.	CHAIR	POST for CUSPIDOR	
		SMALL	BIG
1	TECNODENT ECO	18	5
2	TECNODENT LINDA	18	5
3	JUMBO SDS	18	5
4	PELTON & CRANE	13	5



UNIT PANDA , GALAXY :  
POST FOR LAMP 70 cm

UNIT SUPERTEK :  
POST FOR LAMP 80 cm



# **ETI** DENTAL INDUSTRIES LTD.

New Tel. For Factory and show Room 972-3-5799760, 972-3-5799759  
16, HAKISHON Bney-Brak 51203 ISRAEL  
Mail Address P.O.BOX: 1400 Bney-Brak 51114 Israel.  
FAX. No. 972-3-5799758

## HEIGHT OF POST FOR UNIT & CUSPIDOR FLOOR MOUNT

